

LONDON  
SCHOOL of  
HYGIENE  
& TROPICAL  
MEDICINE



Hillemann, D; Haasis, C; Andres, S; Behn, T; Kranzer, K (2018)  
Validation of the FluoroType MTBDR assay for the detection of ri-  
fampicin and isoniazid resistance in Mycobacterium tuberculosis com-  
plex isolates. *Journal of clinical microbiology*. ISSN 0095-1137 DOI:  
<https://doi.org/10.1128/JCM.00072-18>

Downloaded from: <http://researchonline.lshtm.ac.uk/4647223/>

DOI: [10.1128/JCM.00072-18](https://doi.org/10.1128/JCM.00072-18)

#### Usage Guidelines

Please refer to usage guidelines at <http://researchonline.lshtm.ac.uk/policies.html> or alternatively contact [researchonline@lshtm.ac.uk](mailto:researchonline@lshtm.ac.uk).

Available under license: <http://creativecommons.org/licenses/by/2.5/>

1 **Validation of the FluoroType® MTBDR assay for the detection of rifampicin and**  
2 **isoniazid resistance in *Mycobacterium tuberculosis* complex isolates**

3

4 Doris Hillemann<sup>1</sup>, Carsten Haasis<sup>1</sup>, Sönke Andres<sup>1</sup>, Tobias Behn<sup>1</sup>, Katharina Kranzer<sup>1, 2\*</sup>

5

6 <sup>1</sup>National Reference Laboratory, Research Centre Borstel, Germany

7 <sup>2</sup>London School of Hygiene and Tropical Medicine, United Kingdom

8

9

10

11 Running title: Fluorotype® MTBDR for first line resistance

12 Key words: tuberculosis, molecular methods, isoniazid, rifampicin, drug resistance

13

14

15

16 \*Corresponding author: Katharina Kranzer

17 Mailing address: Forschungszentrum Borstel, Nationales Referenzzentrum für Mykobakterien,

18 Parkallee 18, D-23845 Borstel, Germany

19 Phone: (49)-4537-1887610, Fax: (49)-4537-1883110

20 e-mail: katharina.kranzer@lshtm.ac.uk

21

22

23

24

25

26

27

28

29

30

31

32

33

34

35  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60  
61  
62  
63  
64

## ABSTRACT

For *Mycobacterial tuberculosis* complex the rapid and accurate diagnosis of drug resistance is crucial to ensure early initiation of appropriate therapy. Recently a new molecular diagnostic test, the FluoroType® MTBDR aimed at detecting rifampicin and isoniazid resistance has become available. This study aimed to evaluate the FluoroType® MTBDR in comparison to phenotypic drug susceptibility testing (DST) using *M. tuberculosis* complex isolates. MTBC isolates underwent phenotypic DST and were tested using the FluoroType® MTBDR and Genotype MTBDRplus. Sanger sequencing of the key regions of *rpoB*, *katG*, *inhA* and *aphC* was performed for isolates with discordant phenotypic and molecular results. Furthermore isolates with specific wildtype bands missing in the Genotype MTBDRplus indicating the presence of a mutation were investigated by Sanger sequencing. Specificity and sensitivity defined as the proportion of isolates correctly determined as susceptible and resistant by the FluoroType® TBMDR compared to phenotypic DST were calculated. A total of 180 culture isolates were included; phenotypic DST showed 85 isolates susceptible to isoniazid and rifampicin, 7 with isoniazid mono-resistance, 7 with rifampicin mono-resistance and 81 with multi-drug resistance. Specificity of the FluoroType® MTBDR was 100% (95%CI 96.0-100%) for both rifampicin and isoniazid. Sensitivity was 91.7% (95%CI 83.6-96.6%) for isoniazid and 98.9% (95%CI 93.8-100.0) for rifampicin. The FluoroType® MTBDR has a high sensitivity and specificity for the detection of rifampicin and isoniazid resistance when using culture isolates.

65 **Introduction**

66 Rapid and accurate laboratory diagnosis of *Mycobacterium tuberculosis* drug susceptibility and  
67 resistance are crucial to ensure early initiation of appropriate therapy, to adequately manage  
68 disease and to control further transmission. Conventional drug susceptibility testing (DST)  
69 relies on culture-based methods with results only available after several weeks. Much shorter  
70 turn-around times can be achieved with molecular diagnostics either performed directly on  
71 sputum samples or positive cultures.<sup>1-3</sup>

72 The currently available molecular methods include Xpert® MTB/RIF, line-probe assays  
73 (LPAs), Sanger target sequencing and next generation sequencing.<sup>4-7</sup> LPAs are recommended  
74 by the World Health Organization (WHO) as rapid diagnostic tests for detection of drug  
75 resistance.<sup>8</sup> The technology combines PCR with subsequent reverse hybridization. The  
76 specifically bound amplicons are made visible in a colorimetric detection reaction resulting in  
77 banding patterns suggestive of the presence or absence of mutations in the target genes. The  
78 most widely used LPA, the GenoType MTBDR*plus*, detects the most prevalent mutations in  
79 the *rpoB*, *katG*, *inhA* genes conferring resistance to rifampicin (RMP) and isoniazid (INH). The  
80 GenoType® MTBDR*plus* is designed to specifically detect four *rpoB* mutations (in codons 516,  
81 526 and 531), two *katG* mutations (in codon 315) and four *inhA* mutations in the regulatory  
82 region. Other mutations within the amplified region of the target genes are indicated by the  
83 absence of a wildtype band without the simultaneous presence of a mutation band.

84 Recently a new assay, the FluoroType® MTBDR directed at the same target genes, but using a  
85 different technology has been developed. This assay combines a Linear-After-The-Exponential  
86 (LATE)-PCR<sup>9</sup> together with special probes using lights-on/lights-off detection technology.<sup>10</sup>

87 The readouts of the FluoroType® MTBDR are melting curves. The shapes correspond to  
88 wildtypes or the presence of specific mutations. The FluoroType® MTBDR is an open system,  
89 that identifies characterized mutation *via* a learning software interpreting the melting curves.

90 Practically the FluoroType® MTBDR has several advantages over the GenoType®  
91 MTBDR*plus* i) less hands-on time ii) more rapid results ii) decreased risk of DNA  
92 contamination and iv) automatic interpretation with the possibility to import results directly  
93 into a laboratory information system. However, sensitivity and specificity for detecting RMP  
94 and INH resistance of this new method have not been investigated under routine conditions.

95 This study aimed to compare the results of phenotypic DST, GenoType® MTBDR*plus* and  
96 FluoroType® MTBDR using cultures positive for *M. tuberculosis* complex.

97

98

99 **Methods**

100 Mycobacterial cultures were performed at the National Mycobacterial Reference Laboratory in  
101 Borstel, Germany between May 2016 and October 2017 using mycobacterial growth indicator  
102 tubes (MGIT) (Becton Dickinson, USA) and Loewenstein-Jensen and Stonebrink culture slants  
103 (Artelt-Enclit, Germany). All cultures were identified to species level using the GenoType  
104 MTBC ver 1.X (Hain Lifescience, Nehren, Germany) according to the manufacturer's  
105 instructions. The MGIT 960-Isoniazid-Rifampicin-Ethambutol (IRE) kit (Becton Dickison)  
106 was used according to the manufacturer's instructions. The critical concentrations of  
107 Rifampicin and Isoniazid were 1 µg/ml and 0.1 µg/ml, respectively.<sup>11</sup>  
108 For DNA-extraction an aliquot of 500 µl was obtained from the bottom of a positive MGIT or  
109 two loops (10 µl) of bacterial growth were collected from Loewenstein-Jensen or Stonebrink  
110 culture slants and added to 500 µl of 0.9% NaCl solution. Samples were centrifuged at 3000 x  
111 g for 15 minutes. Supernatants were removed, and DNA extracted using the Fluorolyse® kit  
112 (Hain Lifescience, Nehren, Germany) according to the manufacturer's instructions. In brief,  
113 100 µl of lysis buffer and 2 µl of internal control were added to the remaining pellets. The  
114 mixed samples were incubated at 95°C for 5 minutes, centrifuged for 1 min and resuspended in  
115 100 µl neutralization buffer. Following centrifugation at 10000 x g for 5 minutes supernatants  
116 were transferred into a 0.5 ml cryovial. DNA was stored at -20°C. For each batch, one negative  
117 control was included using 100 µl of lysis buffer, 100 µl of neutralization buffer and 2 µl of  
118 internal control.  
119 DNA was analysed with the Genotype MTBDR*plus* version 2 and FluoroType® MTBDR  
120 assays (Hain Lifescience, Nehren, Germany) as per manufacturer's instructions. In brief, for  
121 the FluoroType® MTBDR 20 µl of DNA-solution was added to 6 µl master mix component A  
122 and 14 µl master mix component B and transferred into 96-well plates. Plates were then sealed,  
123 centrifuged and analyzed via FluoroCycler® 96. For each plate, one positive and one negative  
124 control were included. For the Genotype MTBDR*plus* version 2, 5 µl of DNA-solution were  
125 added to 45 µl of master mix, containing 10µl master mix component A and 35 µl component  
126 B. 20 µl of PCR-products were added to 20µl denaturation reagent in a hybridization well.  
127 Hybridization reagents were prepared, preheated and loaded onto the GT-Blotter® 48 (Hain  
128 Lifescience, Nehren, Germany). After an initial hybridization step, strips were added and  
129 hybridization performed as per standard protocol. For each batch, one negative control (45 µl  
130 master mix, 5 µl DNA free water) was tested. Interpretation of the results was performed  
131 according to the manufacturer's instructions.

132 The package insert of the FluoroType® MTBDR identifies the following *rpoB*, *katG* and *inhA*  
133 mutations: T508A, S509T, E510H, L511P, S512K, Q513L, Q513P, Q513R, D516A, D516F,  
134 D516V, D516Y, N518I, S522L, S522Q, H526C, H526D, H526G, H526L, H526N, H526P,  
135 H526Q, H526R, H526S, H526Y, R529K, S531F, S531L, S531L, S531Q, S531W, L533E,  
136 L533P, S315T1, S315T2, S315N, S315R, G-17T, A-16G, C-15T, G-9A, T-8A, T-8C and T-  
137 8G. In contrast the Genotype MTBDRplus is able to detect and specify the following  
138 mutations: D516V, H526Y, H526D, S531L, S315T1, S315T2, C-15T, A-16G, T-8C and T-8A.  
139 Molecular results were interpreted without knowledge of phenotypic DST results. Sanger  
140 sequencing in the key regions of *rpoB*, *katG*, *inhA* and *aphC*<sup>12</sup> was performed for isolates with  
141 discordant results either between the two molecular methods or the molecular and phenotypic  
142 methods. If the Genotype MTBDRplus showed a missing wildtype band, but not mutation band  
143 Sanger sequencing in the respective key region was performed to confirm the FluoroType®  
144 MTBDR result. Sanger sequencing was performed using an ABI 3130xl genetic analyzer  
145 (Applied Biosystems) and an ABI BigDye Terminator cycle sequencing kit (version 3.1)  
146 according to the manufacturer's instructions.

147 Molecular results were coded as wildtype, individual mutation, mutations in a region of the  
148 gene and indeterminate. Interpretation of molecular results was susceptible, resistant and  
149 indeterminate. Indeterminate results in the Genotype MTBDRplus were defined as detection of  
150 MTBC DNA without the presence of a gene locus controls. All samples yielding an  
151 indeterminate result were re-tested using a new aliquot of DNA. A mixed culture  
152 (heteroresistance) was defined as both wild-type and mutant DNA being present either  
153 identified by wild-type and mutant band (Genotype MTBDRplus) and/or respective double  
154 peaks in the DNA sequence.

155 To determine the accuracy of the FluoroType® MTBDR assay, results were compared with  
156 phenotypic DST results. Sensitivity and specificity were calculated for RMP and INH  
157 separately. Sensitivity was defined as the proportion of isolates correctly determined as  
158 resistant by the FluoroType® MTBDR compared to the phenotypic DST results. Specificity  
159 was defined as the proportion of isolates correctly determined as susceptible by the  
160 FluoroType® MTBDR compared to the phenotypic DST results. McNemar test on paired  
161 proportions was used to compare sensitivity of the FluoroType® MTBDR and the Genotype  
162 MTBDRplus for INH. All statistical analysis were performed using Stata 14.0.

163 Information was extracted from the laboratory information system anonymously. The ethics  
164 committee of the University of Luebeck approved the study. Individual consent was not sought  
165 as no additional patient information was obtained.

166

167 **Results**

168 A total of 180 cultures were included in the study. The majority was identified as *M.*  
169 *tuberculosis* (n=173), followed by *M. africanum* (n=5), *M. bovis ssp. bovis* (n=1) and *M. bovis*  
170 *var. BCG* (n=1). Analysis of RMP and INH resistance included all isolates regardless of  
171 species. 85 isolates tested phenotypically susceptible for both RMP and INH, 7 were INH  
172 mono-resistant, 7 RMP mono-resistant and 81 resistant to both RMP and INH (table 1). The  
173 number of indeterminate results by FluoroType® MTBDR was 1 (0.6%) for *rpoB*, 2 (1.1%) for  
174 *inhA* and 6 for *katG* (3.3%) initially. Following repeated DNA extraction the number of  
175 indeterminate results was 0 (0.0%) for *rpoB*, 1 (0.6%) for *inhA* and 3 for *katG* (1.7%).  
176 Sensitivity of the FluoroType® MTBDR for INH was 91.7% (77/84, 95%CI 83.6-96.6) and  
177 specificity 100% (92/92, 96.1-100.0). For RMP the sensitivity was 98.9% (87/88, 95%CI 93.8-  
178 100.0) and specificity 100.0% (90/90, 95%CI 96.0-100.0). The most frequent *rpoB* mutation  
179 detected by FluoroType® MTBDR was S531L (n=56). A total of 7 *rpoB* mutations remained  
180 unidentified. Mutations were identified by Sanger *rpoB* sequencing as H526L (n=3), S531L  
181 (n=1), Q513K (n=1), D516Y (n=1) and one mixed culture with a H526D and a S531L  
182 mutation. A total of four mixed culture were identified in the study. One of those cultures was a  
183 mixture of an *rpoB* wildtype isolate and an isolate with a L533P mutation. The FluoroType®  
184 MTBDR detected the L533P mutations while the Genotype MTBDR*plus* did not. The  
185 FluoroType® MTBDR identified three isolates as having H526P mutations while Sanger  
186 sequencing revealed H526L mutations. The FluoroType® MTBDR falsely detected an H526R  
187 mutation in a mixed culture with one isolate harboring a H526D and the other a S531L  
188 mutation. Of the seven cultures falsely tested as INH susceptible by the FluoroType® MTBDR  
189 four had *katG* and/or *inhA* mutations detected by the Genotype MTBDR*plus*. This resulted in a  
190 higher, but non-significant sensitivity of the MTBDR*plus* (96.0%; 95%CI 89.9-99.3) for INH  
191 compared to the FluoroType® MTBDR (91.7%; 95%CI 83.6-96.6) (p=0.12).

192

193 **Discussion**

194 This study shows high sensitivity and specificity of the FluoroType® MTBDR for the  
195 detection of RMP and INH resistance using *M. tuberculosis* complex isolates. The proportion  
196 of indeterminate results was small. Repeat DNA extraction resolved 5 of 9 indeterminate  
197 results suggesting that insufficient or low quality DNA might have been the problem. Two of  
198 the isolates that tested repeatedly indeterminate were mixed cultures which might explain the  
199 result.

200 For RMP, the sensitivity and specificity of FluoroType® MTBDR were comparable to those  
201 reported for established molecular diagnostics such as the Genotype MTBDR*plus* and the  
202 Xpert MTB/RIF.<sup>2,3,13</sup> While the presence or absence of an *rpoB* mutation was correctly  
203 detected in all isolates, the FluoroType® MTBDR failed to identify the exact mutation in 10%  
204 (9/87) of cases. Given that silent *rpoB* mutations are rare,<sup>14-16</sup> one might argue that detecting  
205 the presence of an *rpoB* mutation itself is sufficient. In fact, the most widely used molecular  
206 diagnostic test, the Xpert MTB/RIF only indicates the presence of RMP resistance without  
207 information about the associated *rpoB* mutation.

208 The sensitivity of the FluoroType® MTBDR for detecting *katG* and/or *inhA* mutations was  
209 lower compared to the Genotype MTBDR*plus*. Two of the isolates falsely reported as INH  
210 susceptible by the FluoroType® MTBDR had a rare combination of mutations, a S315T1 *katG*  
211 mutation and a C-17T mutation in the *inhA* promotor region, one isolate had a T-8A mutation  
212 in the *inhA* promotor region. Theoretically the FluoroType® MTBDR has the potential for  
213 continuous improvement when new mutations and their respective melting curve properties are  
214 detected and fed back into the analytic software. Whether and how frequently improvement  
215 through software updates will occur remains to be seen.

216 WHO recommends, that commercial molecular LPAs such as the Genotype MTBDR*plus* may  
217 be used for positive sputum smear specimen or a cultured isolate of *M. tuberculosis* complex to  
218 detect resistance to RMP and INH.<sup>8</sup> The Genotype MTBDR*plus* offers benefit compared to  
219 phenotypic DST when used on cultured isolates. Firstly, results are available within 48 hours of  
220 the culture becoming positive in contrast to the 7-14 days to with culture-based DST.<sup>17</sup>

221 Secondly, these assays are not affected by the presence of contamination with bacteria, fungi or  
222 non-tuberculosis mycobacteria, as are liquid culture-based DST systems. The Genotype  
223 MTBDR*plus* and FluoroType® MTBDR detect mutations in the same three target genes.

224 However, there are some distinct advantages of the FluoroType® MTBDR. Except for DNA  
225 extraction, it provides a fully-automated, closed system with a capacity of testing 95 samples in  
226 one batch. This in turn decreases the risk of DNA contamination, which is less problematic  
227 when using DNA extracted from cultured isolates, but may be a significant problem when  
228 testing primary samples. Full automation also means less hands-on time and automated  
229 interpretation of results reducing the risk of analytic errors. Results can be directly imported  
230 into a laboratory information system thus decreasing transcription errors.

231 This study has several strengths and limitations. It was conducted in a National Tuberculosis  
232 Reference Laboratory with standardized procedures for phenotypic and molecular DST and  
233 long-standing expertise with LPA and Sanger sequencing. Interpretation of molecular results



234 was performed blinded to the phenotypic DST results. Half of the isolates included in the study  
235 were resistant to both RMP and INH. However, the majority of isolates harbored the common  
236 S531L *rpoB* and S315T1 *katG* mutations, possibly skewing the sensitivity estimate upwards.  
237 In summary, the FluoroType® MTBDR assay testing culture isolates of *M tuberculosis*  
238 complex has a high sensitivity for detecting RMP resistance, but a lesser sensitivity for  
239 detection of INH resistance. It offers an alternative to the currently WHO-endorsed Genotype  
240 MTBDR*plus* especially for high throughput laboratories. However, validation studies from  
241 other settings are needed to ensure that isolates with a variety of mutations and combinations of  
242 mutations are investigated. Furthermore studies conducted in low- and middle-income  
243 countries with high burden of MDR-TB should be performed to assess the performance of the  
244 new test in those settings.

245

#### 246 **Acknowledgments**

247 We thank K. Ott, B. Schlüter, and A. Witt, Forschungszentrum Borstel, for their excellent  
248 technical assistance; Hain Lifescience GmbH provided two FluoroType® MTBDR and  
249 Genotype MTBDR*plus* assay kits for free. The project was partly funded by the German Centre  
250 for Infectious Research (DZIF). The funder and Hain Lifescience had no role in the design of  
251 the study and collection, analysis, and interpretation of data and in writing the manuscript.

252

## REFERENCES

- 253 1. World Health Organization The use of molecular line probe assay for the detection of resistance to  
254 isoniazid and rifampicin. Geneva, 2016. <http://apps.who.int/iris/handle/10665/250586> last accessed 10.02.2017
- 255 2. Nathavitharana RR, Cudahy PG, Schumacher SG, Steingart KR, Pai M, Denkinger CM. Accuracy of line  
256 probe assays for the diagnosis of pulmonary and multidrug-resistant tuberculosis: a systematic review and meta-  
257 analysis. *The European respiratory journal* 2017; **49**(1).
- 258 3. Steingart KR, Schiller I, Horne DJ, Pai M, Boehme CC, Dendukuri N. Xpert(R) MTB/RIF assay for  
259 pulmonary tuberculosis and rifampicin resistance in adults. *The Cochrane database of systematic reviews* 2014; **1**:  
260 Cd009593.
- 261 4. Theron G, Peter J, Richardson M, Warren R, Dheda K, Steingart KR. GenoType(R) MTBDRsl assay for  
262 resistance to second-line anti-tuberculosis drugs. *The Cochrane database of systematic reviews* 2016; **9**:  
263 Cd010705.
- 264 5. Pankhurst LJ, Del Ojo Elias C, Votintseva AA, Walker TM, Cole K, Davies J, Fermont JM, Gascoyne-  
265 Binzi DM, Kohl TA, Kong C, Lemaitre N, Niemann S, Paul J, Rogers TR, Roycroft E, Smith EG, Supply P, Tang  
266 P, Wilcox MH, Wordsworth S, Wyllie D, Xu L, Crook DW. Rapid, comprehensive, and affordable mycobacterial  
267 diagnosis with whole-genome sequencing: a prospective study. *The Lancet Respiratory medicine* 2016; **4**(1): 49-  
268 58.
- 269 6. Walker TM, Kohl TA, Omar SV, Hedge J, Del Ojo Elias C, Bradley P, Iqbal Z, Feuerriegel S, Niehaus  
270 KE, Wilson DJ, Clifton DA, Kapatai G, Ip CL, Bowden R, Drobniowski FA, Allix-Beguec C, Gaudin C, Parkhill  
271 J, Diel R, Supply P, Crook DW, Smith EG, Walker AS, Ismail N, Niemann S, Peto TE. Whole-genome  
272 sequencing for prediction of Mycobacterium tuberculosis drug susceptibility and resistance: a retrospective cohort  
273 study. *The Lancet Infectious diseases* 2015; **15**(10): 1193-202.
- 274 7. Zignol M, Dean AS, Alikhanova N, Andres S, Cabibbe AM, Cirillo DM, Dadu A, Dreyer A, Driesen M,  
275 Gilpin C, Hasan R, Hasan Z, Hoffner S, Husain A, Hussain A, Ismail N, Kamal M, Mansjo M, Mvusi L, Niemann  
276 S, Omar SV, Qadeer E, Rigouts L, Ruesch-Gerdes S, Schito M, Seyfaddinova M, Skrahina A, Tahseen S, Wells  
277 WA, Mukadi YD, Kimerling M, Floyd K, Weyer K, Raviglione MC. Population-based resistance of  
278 Mycobacterium tuberculosis isolates to pyrazinamide and fluoroquinolones: results from a multicountry  
279 surveillance project. *Lancet Infect Dis* 2016; **16**(10): 1185-92. doi: 10.016/S473-3099(16)30190-6. Epub 2016 Jul  
280 7.
- 281 8. World Health Organisation The use of molecular line probe assays for the detection of resistance to  
282 second-line anti-tuberculosis drugs. Geneva, 2016. <http://www.who.int/tb/WHOPolicyStatementSLLPA.pdf> last  
283 accessed 19.10.2016
- 284 9. Sanchez JA, Pierce KE, Rice JE, Wangh LJ. Linear-after-the-exponential (LATE)-PCR: an advanced  
285 method of asymmetric PCR and its uses in quantitative real-time analysis. *Proc Natl Acad Sci U S A* 2004; **101**(7):  
286 1933-8. doi: 10.073/pnas.0305476101. Epub 2004 Feb 9.
- 287 10. Rice LM, Reis AH, Jr., Wangh LJ. Virtual Barcoding using LATE-PCR and Lights-On/Lights-Off  
288 probes: identification of nematode species in a closed-tube reaction. *Mitochondrial DNA A DNA Mapp Seq Anal*  
289 2016; **27**(2): 1358-63. doi: 10.3109/19401736.2014.947581. Epub 2014 Aug 11.
- 290 11. Cambau E, Viveiros M, Machado D, Raskine L, Ritter C, Tortoli E, Matthys V, Hoffner S, Richter E,  
291 Perez Del Molino ML, Cirillo DM, van Soolingen D, Bottger EC. Revisiting susceptibility testing in MDR-TB by  
292 a standardized quantitative phenotypic assessment in a European multicentre study. *The Journal of antimicrobial*  
293 *chemotherapy* 2015; **70**(3): 686-96.
- 294 12. Sajduda A, Brzostek A, Poplawska M, Augustynowicz-Kopec E, Zwolska Z, Niemann S, Dziadek J,  
295 Hillemann D. Molecular characterization of rifampin- and isoniazid-resistant Mycobacterium tuberculosis strains  
296 isolated in Poland. *J Clin Microbiol* 2004; **42**(6): 2425-31.
- 297 13. Bai Y, Wang Y, Shao C, Hao Y, Jin Y. GenoType MTBDRplus Assay for Rapid Detection of Multidrug  
298 Resistance in Mycobacterium tuberculosis: A Meta-Analysis. *PloS one* 2016; **11**(3): e0150321.
- 299 14. Mathys V, van de Vyvere M, de Droogh E, Soetaert K, Groenen G. False-positive rifampicin resistance  
300 on Xpert(R) MTB/RIF caused by a silent mutation in the rpoB gene. *The international journal of tuberculosis and*  
301 *lung disease : the official journal of the International Union against Tuberculosis and Lung Disease* 2014; **18**(10):  
302 1255-7.
- 303 15. Mokaddas E, Ahmad S, Eldeen HS, Al-Mutairi N. Discordance between Xpert MTB/RIF assay and  
304 Bactec MGIT 960 Culture System for detection of rifampin-resistant Mycobacterium tuberculosis isolates in a  
305 country with a low tuberculosis (TB) incidence. *Journal of clinical microbiology* 2015; **53**(4): 1351-4.
- 306 16. Ocheretina O, Byrt E, Mabou MM, Royal-Mardi G, Merveille YM, Rouzier V, Fitzgerald DW, Pape JW.  
307 False-positive rifampin resistant results with Xpert MTB/RIF version 4 assay in clinical samples with a low  
308 bacterial load. *Diagn Microbiol Infect Dis* 2016; **85**(1): 53-5. doi: 10.1016/j.diagmicrobio.2016.01.009. Epub Jan  
309 15.
- 310 17. Simons SO, van Soolingen D. Drug susceptibility testing for optimizing tuberculosis treatment. *Curr*  
311 *Pharm Des* 2011; **17**(27): 2863-74.
- 312

**Table 1: Results of the FluoroType® MTBDR, GenoType MTBDRplus and phenotypic DST including all 180 isolates**

N	FluoroType® MTBDR					GenoType MTBDRplus					Phenotypic DST		Sanger Sequencing			
	<i>rpoB</i>	<i>katG</i>	<i>inhA</i>	Molecular DST Interpretation		<i>rpoB</i>	<i>katG</i>	<i>inhA</i>	Molecular DST Interpretation		RMP	INH	<i>rpoB</i>	<i>katG</i>	<i>inhA</i>	<i>alpC</i>
				RMP	INH				RMP	INH						
85 <sup>b</sup>	WT	WT	WT	S	S	WT	WT	WT	S	S	S	S				
6	WT	S315T1	WT	S	R	WT	S315T1	WT	S	R	S	R				
1	WT	WT	C-15T	S	R	WT	WT	C-15T	S	R	S	R				
1	WT	S315T1	C-15T	S	R	WT	S315T1	C-15T	S	R	R	R	Mutation outside the hotspot region			
3	S531L	WT	WT	R	S	S531L	WT	WT	R	S	R	S				
1	S531L	WT	WT	R	S	S531L	WT	WT	R	S	R	R	WT	WT	WT	
3	S531L	WT	WT	R	S	530-533	WT	WT	R	S	R	S	S531L			
1	S531L	WT	WT	R	S	530-533	WT	WT	R	S	R	R	S531L	WT	WT	C-10T
1	S531L	WT	WT	R	S	S531L	WT	T-8A	R	R	R	R				T-8A
1	S531L	IND	WT	R	IND	S531L	IND	WT	R	IND	R	R		WT	WT	WT
1	S531L	MUT	WT	R	R	S531L	S315T1	WT	R	R	R	R	S315T1			
2	S531L	S315T2	WT	R	R	S531L	S315T2	WT	R	R	R	R				
2	S531L	S315T1	C-15T	R	R	S531L	S315T1	C-15T	R	R	R	R				
1	S531L	S315T1	C-15T	R	R	530-533	S315T1	C-15T	R	R	R	R	S531L			
1	S531L	S315T1	T-8C	R	R	S531L	S315T1	T-8C	R	R	R	R				
38	S531L	S315T1	WT	R	R	S531L	S315T1	WT	R	R	R	R				
1	S531L	S315T1	WT	R	R	530-533	S315T1	WT	R	R	R	R	S531L			
5	D516V	S315T1	WT	R	R	D516V	S315T1	WT	R	R	R	R				
1	D516V	S315T1	WT	R	R	513-519	S315T1	WT	R	R	R	R	D516V			
1 <sup>a</sup>	D516Y	WT	IND	R	IND	513-519 S531L	IND	C-15T	R	R	R	R	D516Y S531L	WT	C-15T	
1 <sup>c</sup>	D516Y	WT	WT	R	S	513-519	S315T1	-15 to -17	R	R	R	R	D516Y	S315T1	C-17T	

3 <sup>c</sup>	H526D	S315T1	WT	R	R	H526D	S315T1	WT	R	R	R	R				
1	H526D	WT	C-15T	R	R	H526D	WT	C-15T	R	R	R	R				
1	H526D	WT	WT	R	S	H526D	WT	WT	R	S	R	S				
1	H526N	S315T1	WT	R	R	526-529	S315T1	WT	R	R	R	R	H526N			
3	H526P	S315T1	C-15T	R	R	526-529	S315T1	C-15T	R	R	R	R	H526L			
1 <sup>a</sup>	H526R	IND	WT	R	IND	H526D S531L	S315T1	WT	R	R	R	R	H526D S531L	S315T1		
1	H526Y	IND	WT	R	IND	H526Y	S315T1	WT	R	R	R	R	S315T1			
2 <sup>c</sup>	H526Y	S315T1	WT	R	R	H526Y	S315T1	WT	R	R	R	R				
1	H526Y	WT	WT	R	S	H526Y	S315T1	WT	R	R	R	R	S315T1			
1	L533P	S315T1	T-8C	R	R	530-533	S315T1	T-8C	R		R	R	L533P			
1 <sup>a</sup>	L533P	S315T1	WT	R	R	WT	S315T1	WT	S	R	R	R	WT L533P			
3 <sup>c</sup>	MUT	S315T1	WT	R	R	526-529	S315T1	WT	R	R	R	R	H526L			
1	MUT	S315T1	WT	R	R	530-533	S315T1	WT	R	R	R	R	S531L			
1	MUT	S315T1	WT	R	R	510-517	S315T1	WT	R	R	R	R	Q513K			
1 <sup>a</sup>	MUT	WT	WT	R	S	H526D S531L	WT	WT	R	S	R	R	H526D S531L	WT	WT	WT
1 <sup>c</sup>	MUT	WT	WT	R	S	513-519	S315T1	-15 to -17	R	R	R	R	D516Y	S315T1	C-17T	

314

315 DST, drug susceptibility testing; RMP, rifampicin; INH, isoniazid; WT, wildtype; S, susceptible; R, resistant; IND, indeterminate; MUT, unidentified mutation;

316 when wildtype bands were missing in the Genotype MTBDRplus assay the base-pair region is indicated (i.e. 526-529) ;

317 <sup>a</sup>mixed cultures; <sup>b</sup>one *M. bovis* ssp. *bovis* and one *M. bovis* var. BCG included, <sup>c</sup>one *M africanum* included,

318

319

320

321

322